

# **Corrigendum: Augmented COlorimetric NANoplasmonic (CONAN) Method for Grading Purity and Determine Concentration of EV Microliter Volume Solutions**

Andrea Zendrini<sup>1,2\*</sup>, Lucia Paolini<sup>3,4\*</sup>, Sara Busatto<sup>4,5</sup>, Annalisa Radeghieri<sup>3,4</sup>, Miriam Romano<sup>3,4</sup>, Marca H. M. Wauben<sup>6</sup>, Martijn J. C. van Herwijnen<sup>6</sup>, Peter Nejsum<sup>7</sup>, Anne Borup<sup>7</sup>, Andrea Ridolfi<sup>4,8,9</sup>, Costanza Montis<sup>4,9</sup> and Paolo Bergese<sup>2,3,4</sup>

<sup>1</sup> Department of Animal Science, Food and Nutrition-Università Cattolica del Sacro Cuore, Piacenza, Italy, <sup>2</sup> Consorzio Interuniversitario Nazionale per la Scienza e la Tecnologia dei Materiali, Florence, Italy, <sup>3</sup> Department of Molecular and Translational Medicine, University of Brescia, Brescia, Italy, <sup>4</sup> Consorzio Sistemi a Grande Interfase, Department of Chemistry, University of Florence, Sesto Fiorentino, Italy, <sup>5</sup> Department of Transplantation, Mayo Clinic, Jacksonville, FL, United States, <sup>6</sup> Department of Biochemistry and Cell Biology, Faculty of Veterinary Medicine, Utrecht University, Utrecht, Netherlands, <sup>7</sup> Department of Clinical Medicine, Faculty of Health, Aarhus University, Aarhus, Denmark, <sup>8</sup> Istituto per lo Studio dei Materiali Nanostrutturati (CNR-ISMN), Bologna, Italy, <sup>9</sup> Department of Chemistry, University of Florence, Sesto Fiorentino, Italy

# **OPEN ACCESS**

Edited and reviewed by:

lija Kusters, Harvard Medical School, United States

### \*Correspondence:

Andrea Zendrini andrea.zendrini@unibs.it Lucia Paolini lucia.paolini@unibs.it

# Specialty section:

This article was submitted to Synthetic Biology, a section of the journal Frontiers in Bioengineering and Biotechnology

> Received: 01 March 2021 Accepted: 10 March 2021 Published: 13 April 2021

#### Citation:

Zendrini A, Paolini L, Busatto S, Radeghieri A, Romano M, Wauben MHM, van Herwijnen MJC, Nejsum P, Borup A, Ridolfi A, Montis C and Bergese P (2021) Corrigendum: Augmented COlorimetric NANoplasmonic (CONAN) Method for Grading Purity and Determine Concentration of EV Microliter Volume Solutions. Front. Bioeng. Biotechnol. 9:674507. doi: 10.3389/fbioe.2021.674507 Keywords: extracellular vesicles, synthetic vesicles, liposomes, purity, titration, particle number, nanoparticles, nanoplasmonics

# A Corrigendum on

# Augmented COlorimetric NANoplasmonic (CONAN) Method for Grading Purity and Determine Concentration of EV Microliter Volume Solutions

by Zendrini, A., Paolini, L., Busatto, S., Radeghieri, A., Romano, M., Wauben, M. H. M., et al. (2020). Front. Bioeng. Biotechnol. 7:452. doi: 10.3389/fbioe.2019.00452

In the original article, there was an error. The reported concentration of the trisodium citrate dihydrate solution, "3.5 mM," used to synthesize gold nanoparticles (AuNPs), was incorrect. This error inevitably leads to failure in gold nanoparticle synthesis. The correct trisodium citrate dehydrate concentration should be "34.0 mM."

A correction has been made to the **Methods** section, subsection **Preparatory Procedures**, paragraph, under the **AuNP Synthesis** method:

# AuNP Synthesis

AuNPs are synthesized through classic Turkevich's citrate reduction method (Turkevich et al., 1951).

- Dissolve trisodium citrate  $\cdot$  2H<sub>2</sub>O in HPLC grade water to a final concentration of 34.0 mM (1.0% wt).
- Prepare 20 ml of 1 mM HAuCl<sub>4</sub> in HPLC grade water.
- Boil the HAuCl<sub>4</sub> solution under continuous stirring.
- Inject 2 ml of trisodium citrate 2H<sub>2</sub>O 34.0 mM into the boiling solution. **Crucial step:** trisodium citrate injection speed influences AuNP size and therefore final concentration. A one-shot injection leads to smaller particles (usually around 12–15 nm); slow injection will result in NP aggregation and precipitation.
- Keep stirring and wait for the solution to change color from the original pale yellow to wine red.

- After 10 min, cool the solution in a water-ice bath for 5 min.
- Store the AuNPs at 4°C. **Crucial step:** 10–20 nM AuNP solutions are metastable and tend to form aggregates and precipitate with time. In our experience they keep their properties with respect to the assay up to 1 month of storage at the above described conditions. Anyhow it is strongly suggested to freshly (re)determine before any usage the key

# REFERENCES

Turkevich, J., Stevenson, P. C., and Hillier, J. (1951). A study of the nucleation and growth processes in the synthesis of colloidal gold. *Discuss. Faraday Soc.* 11, 55–75. doi: 10.1039/df9511100055 characteristics of the AuNP solution which are relevant for the assay as described in section Determination of AuNP Molar Concentration and Spectral Properties.

The authors apologize for this error and state that this does not change the scientific conclusions of the article in any way. The original article has been updated.

Copyright © 2021 Zendrini, Paolini, Busatto, Radeghieri, Romano, Wauben, van Herwijnen, Nejsum, Borup, Ridolfi, Montis and Bergese. This is an open-access article distributed under the terms of the Creative Commons Attribution License (CC BY). The use, distribution or reproduction in other forums is permitted, provided the original author(s) and the copyright owner(s) are credited and that the original publication in this journal is cited, in accordance with accepted academic practice. No use, distribution or reproduction is permitted which does not comply with these terms.