

Prolonged Animal Observation by Use of Digitized Videodisplays

BERRIE M. SPRUIJT AND WILLEM HENDRIK GISPEN

*Institute of Molecular Biology and Rudolf Magnus Institute for Pharmacology
State University of Utrecht, Utrecht, The Netherlands*

Received 22 March 1983

SPRUIJT, B. M. AND W. H. GISPEN. *Prolonged animal observation by use of digitized videodisplays*. PHARMACOL BIOCHEM BEHAV 19(5) 765-769, 1983.—Quantitative assessment of animal behavior requires not only much time and good training, but is always dependent upon the interpretation of the observer. Therefore we employed a commonly-used microcomputer to analyze, without intervention of a human observer, the picture obtained via videocamera, a digitizer and disk-drive. The program, written in BASIC and listed in Appendix I, generates pictures of the monitored animal in the graphics of the computer. These pictures can be analyzed immediately or stored on disk for further analysis.

Computerized animal observation Pattern recognition Microcomputer Videodigitizer

IN BEHAVIORAL pharmacology one often deals with the problem of quantifying certain clearly described behaviors in different groups of animals. For the purpose of screening drugs, several typical and easily recognized behaviors are used such as akinesia, hyperactivity, hypoactivity, stereotypic gnawing, licking, biting, wet dog shake and shivering.

Often special electronic equipment adapted for particular tests has been developed to measure easily quantifiable parameters; which are derived from animal movement [5]. One disadvantage of automation by making the appropriate hardware is the limited flexibility of the equipment. Moreover, behavioral elements—like grooming—which are more defined by attitude than movement, are more difficult to measure. Often the intervention of a human observer remains a prerequisite for registering animal behavior. This intervention not only introduces bias but is also dependent on sustained attention of the observer and might therefore interfere with the scoring consistency.

Here a system is described which can replace a time sampling method for observing two animals at the same time. In our laboratory it is used for observing the grooming behavior of two rats. Both spontaneously occurring grooming as well as peptide-induced grooming are prolonged behavioral patterns [7]. As a consequence demonstrating differences in grooming behavior requires a lot of observation time and suffers from the afore-mentioned disadvantages inherent in a human observer.

METHOD

The microcomputer used is an Apple II Plus with Applesoft BASIC and 48K RAM. Other computers compatible with Apple will do as well. The Computer is equipped with a disk-drive and a graphic printer. An essential interface is the Apple videodigitizer. This interface is connected to a standard videocamera and generates digitized pictures every 4 (or

8) seconds of the display of the camera screen. This picture will be stored in the graphics memory of the computer. The digitizer also contains microswitches for regulating contrast and brightness. The animals used for observation are housed individually in cages of Plexiglas (dimensions: 14.5×14.5×25 cm) with a black background for the required contrast with white rats [3].

Reliability is demonstrated in a small experiment in which a human observer and our proposed system were both analyzing animal behavior. Four times two rats—one experimental and one control—were observed for two hours.

The experimentals were treated with 0.3 µg/3 µl ACTH to induce the display of excessive grooming, controls received 3 µl of saline.

Software

The manual of the digitizer explains how to make a picture and display it on the screen. Different strategies for further processing the pictures can be followed: (1) immediate analysis of each picture or (2) storage on disk and application of a special program for future analysis of these pictures. Of course, storage is preferable, however, the number of pictures is restricted to the disk storage capacity available. When Apple DOS 3.3 system is used 15 pictures (8,192 bytes per picture equals 33 sectors on disk) can be stored on one disk (single sided, single density). Pattern analysis is relatively time consuming, so immediate analysis results in a reduced sampling rate, since analysis time dictates the interval between pictures analyzed.

The program listed in Appendix I is completely written in BASIC and therefore the analysis of the generated picture also requires some time.

However this can be limited by using a subroutine in assembly language for screening every point in graphics. For understanding the principles of the Apple Graphics and the way our program makes use of it, a program written in

TABLE 1

BEHAVIORAL ANALYSIS PERFORMED BY A HUMAN OBSERVER
AND A COMPUTER SYSTEM

Behavior	Human Observer		Computer	
	E-C		E-C	
Grooming	139 ± 5	63 ± 5*	112 ± 18	56 ± 12*
Sitting	39 ± 5	37 ± 7	58 ± 18	53 ± 10
Sleeping	27 ± 8	119 ± 7*	24 ± 6	112 ± 8*
Rearing	71 ± 12	59 ± 8	81 ± 16	56 ± 12

The behavior of two groups of four animals was analyzed both by a human observer and the proposed computer system. **t*-test revealed significance between the two groups ($p < 0.05$). C: saline-treated; E: ACTH-treated.

BASIC seems more appropriate and is therefore listed and used for the results in Table 1.

In analyzing the two-dimensional display two main strategies can be followed. First, one can give certain examples of different behavioral elements and let the computer compare the digitized pictures with the given pictures; the example with the closest resemblance to the monitored behavior determines the classification. A "rest" category can be used if no resemblance with any example is found. The more graphics pages the computer has the more attractive this approach is (an Apple has only two pages for its graphics (8,192 bytes each).

Secondly, one can transform the digitized pictures to easily recognizable geometrical figures deduced from the posture of the behavior in study. The program shown divides the screen into two halves each containing a picture of an animal. From both pictures the height and the width of the animal at half its height are measured. The width of the previous picture remains available as a cue for classifying behavior and is used for determining movement.

According to the results of Traber [8] and Bolles [1] most time is spent in sitting, grooming, sleeping and exploration. It is our purpose to discriminate between these four behavioral categories. Classification of animal's posture is performed according to the height, the width and the width of the fore-going posture. As can be seen in Fig. 1, which demonstrates three pictures from the screen, rearing and sleeping are rather easy to distinguish. Rearing is characterized by a large height and a small width; sleeping is marked by a small height and a large width. The postures with intermediate values for height and width, grooming and sitting, are more difficult to distinguish; both width and height can be similar in the two postures, therefore a parameter for movement is used as a discriminant. A sitting animal has the same width over the different samples, while a grooming animal is constantly moving. When the animal is not rearing or sleeping the width and the height of that scan and that of the preceding scan are used for discriminating grooming and sitting. A boundary range is defined in case the height is close to the set point of rearing or sleeping (see Fig. 2). The program uses the height as a main key; when the animal is obviously sleeping or rearing—respectively height larger or smaller than respectively set point plus or minus the range—the width is not taken into account. This range is variable and can be changed in the program (see Appendix I,




FIG. 1. Three examples of monitored behavior of two animals in two cages. The numbers represent the calculated values for height and width according to the computer. Behavior was determined for both animals as is shown. The black spot on the head of the animal is caused by the dental cement for cannula implantation.




FIG. 2. The lines present the heights which determine classification of animal behavior. If the height falls within a boundary range (dotted line), the width is also taken into account (see text).

variable UN). When the posture indicates grooming or sitting the width of the preceding picture is also used for discrimination. The set point must be given for every run of the program, so one can always adapt the values according to the position of the camera with respect to the cage, which is very crucial in this set up. The first part of the program provides projection of the given set points on the monitor, so you can check and, if necessary, adapt them to your requirement.

The program counts every point which makes up the picture on the screen—or every second or third point depending on the scanning rate—as either white or black.

The scanning rate of the screen is varied by altering the step size; the bigger the step size the faster and more inaccurate the scanning performed.

When the strategy of comparing the picture with standard figures is followed a good device for getting these standard figures into the computer is a graphics tablet. The postures can be drawn by hand on the tablet and can be placed in graphics automatically. (For further information see the manual of the tablet.)

RESULTS

Table 1 demonstrates how the program can distinguish saline-treated animals from ACTH-treated animals. ACTH-

treated animals groomed significantly more than the saline-treated animals. This is validated by the data of the human observer. Although the drug-induced phenomenon is recognized, considerable difference between our proposed system and the human observer appears with respect to the discrimination of sitting and excessive grooming.

DISCUSSION

The results indicate some shortcomings of the proposed program. Rearing and lying down are well discriminated, however a conservative error is made when grooming behavior—performed without any movement—is scored as sitting. When the occurrence of grooming behavior is studied in prolonged behavioral session, the program works sufficiently. In our laboratory it is also used for the quantification of exploratory behavior; in this case grooming and sitting are lumped together. The advantage of this procedure lies not yet in the accuracy but in the possibility of continuous animal observation for instance in combination with continuous drug administration (osmotic micropumps). Moreover

the bias, which always accompanies the human interpretation of animal behavior is incorporated in a program and therefore accessible to discussion. Exchange of programs could be in the onset to more homogeneity in the way animal behavior is observed and interpreted. The equipment described has already often been used for other purposes in psychopharmacology [2, 4, 6]. When test cages are similar and compatible computers are available, minimum of investment is required to take advantage of the listed program. It is realized that the program still is far from perfect and needs modification or improvement if used in other experimental settings. Hopefully it inspires other programmers to add a required subroutine and to discuss the outcome of this type of programs with each other.

ACKNOWLEDGEMENTS

The authors wish to thank Mr. J. Brakkee and Mr. H. v.d. Wildt for their technical assistance, Mr. E. Kluis for the photographs and Mrs. L. Claessens for her help in preparing the manuscript.

REFERENCES

1. Bolles, R. Grooming behavior in the rat. *J Comp Physiol Psychol* **53**: 306–310, 1960.
2. Carroll, M. E., P. A. Santi and L. R. Rudell. A microcomputer system for the control of behavioral experiments. *Pharmacol Biochem Behav* **14**: 415–417, 1980.
3. Gispen, W. H., V. M. Wiegant, H. M. Greven and D. de Wied. The induction of excessive grooming in the rat by intraventricular application of peptides derived from ACTH. Structure-activity studies. *Life Sci* **17**: 645–652, 1975.
4. Hendrie, C. A. and S. Bennett. A microcomputer technique for the detailed analysis of animal behaviour. *Physiol Behav* **30**: 233–235, 1983.
5. Hoglin, A. U., J. E. Hagglund and B. J. Meyerson. A videointerface for behavioral recordings with applications. *Physiol Behav* **30**: 489–492, 1983.
6. Murray, R. B., D. E. Gmerek, A. Cowan and R. J. Tallarida. Use of programmable timer and data logger in the monitoring of animal behavior. *Pharmacol Biochem Behav* **15**: 135–141, 1981.
7. Spruijt, B. and W. H. Gispen. ACTH and grooming behavior in the rat. In: *Hormones and Behaviour in Higher Vertebrates*, edited by R. Gilles. Bielefeld: in press, 1983.
8. Traber, J., H. R. Kleine and W. H. Gispen. Actions of antidepressant and neuroleptic drugs on ACTH- and novelty-induced behavior in the rat. *Eur J Pharmacol* **80**: 407–414, 1982.

APPENDIX

```

1 HOME : VTAB 1
3 PRINT : PRINT "FOR UNDERSTANDING THE CALCULATIONS OF HEIGHT AND WIDTH READ APPLES REFERENCE MANUAL FOR DETAILS ABOUT HIGHT RESOLU
TION GRAPHICS."
4 PRINT "PROVIDE ENOUGH CONTRAST BETWEEN OBJECT AND BACKGROUND"
5 PRINT : PRINT "UNDER REM STATEMENTS THE MEANING OF USED VARIABLES IS EXPLAINED."
6 PRINT : PRINT "PLACE THE CAMERA AT THE SAME HEIGHT AS THE OBJECTS IN SUCH A WAY THAT THE OBJECTS DON'T CROSS THE MIDLINE OF THE S
CREEN"
7 PRINT : PRINT "FOR FINISHING PRESS RESET"
8 FOR I = 1 TO 15000: NEXT I
10 HOME : VTAB 21
20 D$ = CHR$(4): PRINT D$ "PRE2": PRINT "!*": PRINT D$ "PRE0"
30 HPLOT 0,0 TO 0,159: HPLOT 140,0 TO 140,159: HPLOT 260,0 TO 260,159
40 HPLOT 0,0 TO 260,0: HPLOT 0,159 TO 260,159: FOR I = 0 TO 23: HX = I * 8: HPLOT 2, HX: NEXT I
42 FLASH : PRINT "CALIBRATING CAMERA AND DIGITIZER"
43 NORMAL : PRINT "PICTURE OKE ? (Y/N) " : GET Y$
44 IF Y$ < > "Y" THEN PRINT D$ "PRE0": GOTO 10
45 HOME : VTAB 21
50 INPUT "MAXIMUM WIDTH REARING": R$: HPLOT 10,0 TO 10,10: RP% = 7 * R$: HPLOT 10 + RP%,0 TO 10 + RP%,10
60 INPUT "MAXIMUM WIDTH GROOMING": G$: HPLOT 10,0 TO 10,20: GP% = 7 * G$: HPLOT 10 + GP%,0 TO 10 + GP%,20
70 INPUT "MAXIMUM HEIGHT GROOMING " : P$: HPLOT 2,8 * (24 - P%) TO 5,8 * (24 - P%)
80 INPUT "MAXIMUM HEIGHT SLEEPING " : S$: HPLOT 6,8 * (24 - S%) TO 8,8 * (24 - S%)
90 INPUT "YOU WANT TO GIVE NEW DIMENSIONS ? (Y/N)": Y$: IF Y$ = "Y" THEN GOTO 10
100 UN = 1
105 REM UN=UNCERTAINTY RANGE
110 HGR : D$ = CHR$(4)
120 NZ = 2
125 REM NZ REFERS TO STEP-SEIZE; THE BIGGER NZ, THE FASTER THE PICTURE IS SCANNED AND THE LESS ACCURATELY
126 REM R REFERS TO THE SCREEN HALF
130 PRINT D$ "PRE2"
140 PRINT "!*"
150 PRINT D$ "PRE0"

```



```

160 X = 0:N = 0:R = 0:
170 FOR T1 = 0 TO 2
180 TELLER = 0:T8 = 0
190 FOR T2 = T1 * 40 + 9192 TO 9098 + T1 * 40 STEP 128
200 TELLER = TELLER + 1
210 FOR T3 = 1 + R TO 20 + R STEP N%
220 X = PEEK (T2 + T3)
230 IF X > 0 THEN N = N + 1
240 IF X > 0 AND N > 1 THEN GOSUB 760
245 REM N REFRES TO THE NUMBER OF POINTS WHICH MUST BE DETECTED BEFORE DECIDING THE OBJECT IS FOUND
250 IF X > 0 AND (R + N) > 21 THEN N = 0: GOTO 560
260 IF X > 0 AND N > 1 THEN N = 0: GOTO 330
270 NEXT T3
280 NEXT T2
290 NEXT T1
300 IF R = 20 THEN GOTO 560
310 REM DETERMINING HEIGHT LEFT ANIMAL 1000-1010
320 REM CLASSIFYING BEHAVIOR RIGHT ANIMAL
330 RH = T2 + T3 - (9192 + (T1 * 40) + T3)
340 RH = RH / 128:AH = (2 - T1) * 8 + 8 - RH
342 HL = AH
345 REM AH IS THE HEIGHT OF THE LEFT ANIMAL
360 PRINT "HEIGHT LEFT ANIMAL"AH
370 PRINT "WIDTH LEFT ANIMAL NOW AND BEFORE "BL" "LV
380 IF AH < S% - UN THEN PRINT "LEFT ANIMAL IN A HORIZONTAL POSITION":L1 = L1 + 1: GOTO 520
390 IF AH > P% + UN THEN PRINT "LEFT RAT REARS":R1 = R1 + 1: GOTO 520
400 IF (S% - UN) < AH AND AH < S% THEN GOTO 420
410 GOTO 450
420 IF BL > G% THEN PRINT "LEFT ANIMAL IN IN A HORIZONTAL POSITION":L1 = L1 + 1: GOTO 520
430 IF ABS (LV - BL) > 0 THEN PRINT "LEFT ANIMAL GROOMS":G1 = G1 + 1: GOTO 520
440 IF ABS (AH - HL) > 0 THEN PRINT "LEFT ANIMAL IS GROOMING":G1 = G1 + 1: GOTO 520
445 PRINT "LEFT ANIMAL IS SITTING":S1 = S1 + 1: GOTO 520
450 IF AH < P% - UN AND ABS (BL - LV) > 0 THEN PRINT "LEFT ANIMAL GROOMS":G1 = G1 + 1: GOTO 520
460 IF AH < P% - UN AND ABS (AH - HL) > 0 THEN PRINT "LEFT ANIMAL GROOMD":G1 = G1 + 1: GOTO 520
465 IF AH < P% - UN THEN PRINT "LEFT ANIMAL SITS":S1 = S1 + 1: GOTO 520
470 IF AH < P% AND AH > P% - UN THEN GOTO 490
480 GOTO 510
490 IF BL% > R% AND ABS (BL - LV) > 0 THEN PRINT "LEFT ANIMAL GROOMS":G1 = G1 + 1: GOTO 520
500 IF BL% > R% AND ABS (AH - HL) > 0 THEN PRINT "LEFT ANIMAL GROOMS":G1 = G1 + 1: GOTO 520
505 IF BL% > R% THEN PRINT "LEFT ANIMAL SITS":S1 = S1 + 1: GOTO 520
510 PRINT "LEFT RAT REARS ":R1 = R1 + 1
520 TELLER = 0:N = 0:R = 20
530 GOTO 170
540 REM DETERMING HEIGHT RIGHT ANIMAL 2000-2010
550 REM CLASSIFYING BEHAVIOR ACCORDING TO HEIGHT AND WIDTH
560 RH = T2 + T3 - (9192 + (T1 * 40) + T3)
570 RH = RH / 128:AH = (2 - T1) * 8 + 8 - RH
572 HR = AH
575 REM AH REFERS TO THE HEIGHT OF THE RIGHT ANIMAL
580 PRINT "HEIGHT RIGHT ANIMAL"AH
590 PRINT "WIDTH RIGHT ANIMAL NOW AND BEFOPE "BR" "RV
600 IF AH < S% - UN THEN PRINT "RIGHT ANIMAL IN A HORIZONTAL POSITION":L2 = L2 + 1: CALL 62450: GOTO 130
610 IF AH > P% + UN THEN PRINT "RIGHT ANIMAL REARS":R2 = R2 + 1: CALL 62450: GOTO 130
620 IF AH < S% AND AH > S% - UN THEN 640
630 GOTO 670
640 IF BR > G% AND ABS (RV - BR) > 0 THEN PRINT "RIGHT ANIMAL GROOMS":G2 = G2 + 1: CALL 62450: GOTO 130
650 IF BR > G AND ABS (AH - HR) > 0 THEN PRINT "RIGHT RAT GROOMS":G2 = G2 + 1: CALL 62450: GOTO 130
655 IF BR > G% THEN PRINT "RIGHT RAT IS SITTING":S2 = S2 + 1: CALL 62450: GOTO 130
660 PRINT "RIGHT ANIMAL REARS":R2 = R2 + 1: CALL 62450: GOTO 130
670 IF AH < P% - UN AND ABS (RV - BR) > 0 THEN PRINT "RIGHT ANIMAL GROOMS":G2 = G2 + 1: CALL 62450: GOTO 130
680 IF AH < P% - UN AND ABS (AH - HR) > 0 THEN PRINT "RIGHT ANIMALS GROOMS":G2 = G2 + 1:5": CALL 62450: G
OTO 130
685 IF AH < P% - UN THEN PRINT "RIGHT ANIMAL IS SITTING":S2 = S2 + 1: CALL 62450: GOTO 130
690 IF AH < P% AND AH > P% - UN THEN GOTO 710
700 GOTO 730
710 IF BR > R% AND ABS (RV - BR) > 0 THEN PRINT "RIGHT ANIMAL GROOMS":G2 = G2 + 1: CALL 62450: GOTO 130
720 IF BR > R% AND ABS (AH - HR) > 0 THEN PRINT "RIGHT ANIMAL IS GROOMING":G2 = G2 + 1: CALL 62450: GOTO 130
725 IF BR > R% THEN PRINT "RIGHT ANIMAL IS SITTING":S2 = S2 + 1: CALL 62450: GOTO 130
730 PRINT "RIGHT ANIMAL REARS":R2 = R2 + 1
740 PRINT D% "PREO"
750 GOTO 130

```



```
760 IF T1 = 0 THEN T8 = T2 + 40
770 IF T1 = 1 AND TELLER < 5 THEN T8 = 9000: GOTO 790
780 T8 = 9128
785 POKE T8,1
790 N2 = 0
800 IF R = 20 THEN LV = BR: GOTO 820
810 RV = BL
820 REM T8 IS THE HEIGHT AT WHICH THE WIDTH IS DETERMINED
830 FOR T9 = 1 + R TO 20 + R
840 Y = PEEK (T8 + T9)
850 IF Y > 0 THEN N2 = N2 + 1
860 IF Y < > 0 AND N2 < 2 THEN B1 = T9
870 IF Y = 0 AND B1 < > 0 THEN B2 = T9: GOTO 890
880 NEXT T9
890 REM RV AND LV ARE RESPECTIVELY THE WIDTH OF RIGHT AND LEFT ANIMALS
900 BR = B2 - B1:BL = BR:B1 = 0:B2 = 0: RETURN
902 REM STOP BY PRESSING RESET
903 REM TYP GOTO 910 FOR FINAL RESULTS
910 PRINT "REARING LEFT ANIMAL "R1: PRINT "REARING RIGHT ANIMAL "R2: PRINT
911 PRINT "GROOMING LEFT ANIMAL "G1: PRINT "GROOMING RIGHT ANIMAL "G2: PRINT
912 PRINT "SITTING LEFT ANIMAL "S1: PRINT "SITTING RIGHT ANIMAL "S2: PRINT
913 PRINT "SLEEPING LEFT ANIMAL "L1: PRINT "SLEEPING RIGHT ANIMAL "L2
1000 END
```